

personnel in a laboratory environment. For maximum protection against UVA energy, wear UV safety glasses & avoid looking directly at UV light source.



#### Display

| , |      |        |       |
|---|------|--------|-------|
|   | RED  | YELLOW | GREEN |
|   | 30   | 25     | 26    |
|   | BLUE | WHITE  | UV    |
|   | 18   | 32     | 3     |

This analyzer has been designed to determine if tinted lenses meet EC standards for traffic signal safety as outlined in ISO 14889, ISO DIS 8980-3 and ISO 13666.

5 Filter Level LED's are located on the left of the

One LED comes on and stays on. This indicates the filter level that the lens has achieved. (There is no violation of the safety standards).

Fashing LED. As the main display cycles through the colors as they are being analyzed, a LED may change from being steadily lit to flashing. This indicates a safety standard violation. The LED will start to flash as the main display reaches the color that is out of range. As the main display reaches the end of its cycle, the LED stops flashing. It will resume flashing as it reaches the out of range color during subsequent cycles. color during subsequent cycles

NO LED's illuminated. No LED will be lit if the ratio of UV to visible light is too great to meet the standard.

# **Replacement Parts**

Visible Lamp: BPI#60302

1 Amp, 250 Volt Fuse: BPI #59905

#### **Specifications**

The Traffic Signal Safety Spectrum Analyzer™ with RS-232 port (BPI# 211951) is an invaluable aid for quality control of lenses tinted with visible light treatments. The meter may be used for fashion tint control and for critical density adjustments on therapeutic tints such as BPI Diamond Dye™ 500/550 and Diamond Dye<sup>™</sup> 540. It may also be used to check for UVA blocking. It may also be used for traffic signal color density checks, and general filter level checks.

Variations in density and hardness of CR-39<sup>™</sup> lenses typically affect the ability to accept dye. Two lenses that have been in the same dye tank the same amount of time may not come out with the same tint - either in density or color balance.

The meter's digital display indicates the percentage of light passing through a lens in the red (610nm), yellow (580nm), green (558nm), and blue (470nm) range. It also displays a white light (visible) reading as well as a UV reading. It is a quick and accurate way to check the transmission characteristics of lenses.

The system requires 115 volt -220 volt , 50/60 Hz and is fuse protected by a 1 amp, 250 volt glass fuse. Components are UL and CSA recognized.

**NOTE:** Always use the ground wire on the power cord for safe operation; do not ever bypass it.

| HEIGHT         | WIDTH    | LENGTH                     | VOLT     | AGE                                 | WEIGHT           | FUSE                               | AMPERAGE     |
|----------------|----------|----------------------------|----------|-------------------------------------|------------------|------------------------------------|--------------|
| 5 in.          | 6.25 in. | . 6.75 in.                 | 115 or 3 | 220 v.                              | 6 lbs            | 1 amp/250v.                        | 1 amp        |
| 12.7 cm        | 15.87 cn | m 17.14 cm                 |          |                                     | 2.72 kg          |                                    |              |
| LENS CLEARANCE |          | TEST RANGE                 |          | THE SE                              | T-UP KIT INCLUD  | ES THE FOLLOWIN                    | NG PRODUCTS: |
| 0.75 in.       |          | 350 nm to 400 nm (UVA)     |          | Calibration lens     Patient brochu |                  | ochures                            |              |
| 19.05 mm       |          | 400 nm to 700 nm (Visible) |          | • In                                | struction manual | <ul> <li>Patient broken</li> </ul> | ochure stand |

#### Unpacking

When unpacking your instrument, please check to ensure that no concealed damage occurred in transit. If such is noted, save the shipping carton and immediately notify the shipping company's damage control inspector in your area so a claim may be processed. Failure to do this may void any future claim and replacement. Also, call BPI Customer Service so arrangements for a replacement may be made.

# **Setting Up**

To set up your Traffic Signal Safety Spectrum Analyzer™ with RS-232 port, just connect the power cord to a standard outlet convenient to your work area but away from the immediate vicinity of the lens coloring operation. For computer use connect an RS-232 cable from the meter to the serial port on a P.C.

Although BPI's meters are stable and sturdy, they may be adversely affected by excessive humidity and heat.

# Operation

1. The meter has an ON/OFF switch (in back), computer CALIBRATE and CONTINUE buttons. Turn the unit ON. Allow ten minutes after turn on for UV bulb to stabilize.

2. When this unit has stabilized, press CALIBRATE, Place the calibration filter in the light path and press CONTINUE. When instructed to do so, remove the filter and press continue. You should get readings of approximately 100% for all color channels. If not, re-calibrate.

3. Place the lens to be tested over the light aperture (see picture of unit for example), press CONTINUE and the transmission readings will appear. The LED's on the left will indicate the filter level and will flash when a color whose level is unacceptable for traffic signal purposes is detected. Lenses with moderate to high power may cause erroneous readings.

4. Five minutes after calibration the unit will shut itself off; push the CALIBRATE button and you are back to step 2.

5. In the event of any display or computer malfunction, press calibrate to restore normal operation. 6. Use with a Computer: The RS-232 port and the cable are designed to interconnect the BPI Traffic Signal Spectrum Analyzer™ with a P.C. A CD - ROM is supplied with the unit to allow them to communicate.

### **PC System Requirements**

System Requirements: IBM PC-compatible computer. 486 processor or better. 12MB memory. 5 MB disk space. RS-232 serial port. SVGA (800 x 600) display, preferably color. Mouse or other pointing device, CD-ROM drive, (For installation), Windows 95/98, or NT4.0

### Installation of Software:

 Insert the CD-ROM into the CD-ROM drive • From the START menu, select RUN • Enter (CD-ROM drive letter):setup.exe. Choose OK and follow the instructions in the installation program. Re-booting may be necessary to complete the installation. 'READ' and 'Calibrate' may be executed from your PC. All information can be saved and printed out to your PC's printer.

#### Note

The UV reading is 100% of the upper 31% of the 315 to 380nm near UV range. To obtain an approximate value for that range multiply the UV reading by .31.

#### **Filter Levels**

This Analyzer has been designed to determine if tinted lenses meet EC standards for traffic signal safety as outlined in ISO 14889. ISO DIS 8980-3 and ISO 13666. 5 Filter Level LED's are located on the left of the instrument. One LED comes on and stays on. This indicates the filter level that the lens has achieved. (There is no violation of the safety standards). Flashing LED indicates a safety standard violation.

| ISO/DIS 8980-3. Categories for luminous transmittances and the related<br>permissible transmittances in the ultraviolet solar spectral range |   |            |  |  |  |  |  |
|--|---|------------|--|--|--|--|--|
| CATEGORIES   | VISIBLE SPEC                                  | TRAL RANGE | ULTRAVIOLET SPECTRAL RANGE                                 |  |  |  |  |
|  | Range of luminous transmittance ( <i>Tv</i> ) |            | Maximum value of solar UV-A transmittance ( <i>Tsuva</i> ) | Maximum value of solar UV-B transmittance ( <i>Tsuvb</i> ) |  |  |  |
|  | From over                                     | То         | Over 315nm to 380nm UV-A                                   | Over 280nm to 315nm UV-B                                   |  |  |  |
| 0  | 80%   | 100%       | Τν   | Τν   |  |  |  |
| 1  | 43%   | 80%        | Τν   | 0.125 <b>T</b> v   |  |  |  |
| 2  | 18%   | 43%        | Τν   | 0.125 <b>T</b> v   |  |  |  |
| 3  | 8%  | 18%        | 0.5 <b>T</b> v   | 0.125 <b>T</b> v   |  |  |  |
| 4  | 3%  | 8%         | 0.5 <b>T</b> v   | 1.0% absolute  |  |  |  |

|        | ISO 14889. Relative visual att<br>Spectacle lenses conforming to<br>a relative visual atter |
|--------|---|
| RED    | 0.8 (   |
| YELLOW | 0.8 (   |
| GREEN  | 0.6 (   |
| BLUE   | 0.4 (   |
|        |   |

© 2000 BPL All specific product names mentioned herein are trademarks of Brain Power Incorporated, Miami, Florida, USA. (Unless otherwise stated). BPI is a registered trademark with the US Patent Office and v offices in other countries. MANUAL FILE# M2086

NO LED's illuminated. No LED will be lit if the ratio of UV to visible light is too great to meet the standard.

tenuation quotient for signal light recognition ISO 14889 for daylight and night use shall have enuation quotient Q of not less than: (80% of the visible light transmission) (80% of the visible light transmission) (60% of the visible light transmission) (40% of the visible light transmission)

Traffic Signal Spectrum Analyzer with