

BPI® LABORATORY EQUIPMENT

BPI® Vent Hood™

BPI# 15002 (115v.)
BPI# 215002 (220v.)

Keep your area safe and comfortable. Our Vent Hood helps exhaust laboratory fumes, odors and steam with a whisper-quiet turbo fan. The BPI® Turbo Safety Vent Hood™ includes an activated charcoal filter system. Durable construction with easy care finish and stainless steel base. The roomy interior accommodates most BPI® tinting systems, as well as other systems. The unit comes preassembled. It may also be vented to the outside. 3" diameter outlet moves 900 cfm of free air.



HEIGHT	WIDTH	LENGTH	VOLTAGE	WEIGHT	FUSE	AMPERAGE
42 in.	36 in.	28 in.	110v or 220v.	200 lbs.	10 A/ 250 v.	10 amps
107 cm	91.5 cm	71.12 cm		90.71 kg		

BPI® Production Vent Hood™

BPI# 15001 (115v.)
BPI# 215001 (220v.)

Two 3" diameter outlets move 900 cfm of free air each.

HEIGHT	WIDTH	LENGTH	VOLTAGE	WEIGHT	FUSE	AMPERAGE
46 in.	51 in.	28 in.	110v or 220v.	200 lbs.	10 A/ 250 v.	10 amps
116.84 cm	129.5 cm	71.12 cm		90.71 kg		

BPI® Vapor Lip Exhaust System™



Fits many different kinds of tinting systems

BPI® Safety Vapor Exhaust System™ removes fumes from above your BPI® lens tinting instrument. The laminar flow design of the extractor head creates an air curtain above the tinting instrument. The unit features a powerful vacuum pump that can be located up to twelve feet away from the tinting instrument. (The vacuum pump can then be vented to the outside or the air can be re-circulated into the lab after passing through the optional charcoal filters). Available for BPI® Super Color 9™, BPI® Solar Color 6™, BPI® Solar Color 4L™, BPI® Mini Tank 8™.

USA Toll-free: 1-800-327-2250 USA International Phone: 001-305-264-4465 UK Toll-free: 0800 626278 UK International Phone: 01788-568686 France: APPEL GRATUIT 0800 90 81 57 Mexico: LLAMADA GRATUITA: 95 800 225 5274 Italy: NUMERO VERDE: 1678 72037 Germany: RUFEN SIE UNS AN ZUM NULLTARIF: 0800 181-8058 Belgium: NUMERO VERT: 0800 73703